


Archwire Seating Instrument



PART NUMBER

246.1002

Used To Fully Seat The Archwire Into The Archwire Slot.

The Archwire Seating Instrument is the perfect tool for fully seating and de-torque larger cross-section archwires like .019" x .019" and .020" x .020". The tool also simplifies the installation of archwires with smaller cross-sections and more resilient materials like Nickel-Titanium. The ligating door on the Pitts21 or Clear21 brackets are easily closed once the archwire is fully seated into the archwire slot.

Features & Benefits

- Fully seats and de-torques larger cross-section archwires
- Simplifies installation stiffer archwire materials
- Helps for easier ligating door closure